Spray Dryer Pulvis Mini Spray

Spray Dryer (For Granulating, Drying, Mixing)


GB210-B



50g~300g


Temp. contr

40°C~ 220°C



Spray nozzle (selectable)

Tow-way



Display language

Japanese, English, Chinese

Designed to granulate powder and dry wet powder using a fluid bed. This is a fluid bed drying granulator used in combination with the basic unit GB210 and Mini-bed attachment GF200.

- Conditions such as hot air temperature, air amount, binder liquid flow amount can be set easily with the setting dial on the front of the unit
- The chamber is made of ultra hard glass and the user can observe status of the fluid bed or spraying status. Also, the flowage meter, the spraying pressure meter, the chamber inlet/outlet temperature indicator are useful for evaluation of data
- The unit can also be used as a spraying dryer by installing the mini spray attachment GF300 (optional)
- The unit has an automatic lift as a standard to enable convenient installation or removal of the glass chamber attachment

Control Panel



Inlet temperature, outlet temperature, and drying air amount are digitally displayed. Setting is made on the touch panel that allows operation settings, operation status displayas well as error display, and settings of various operation conditions. as well as error display, and settings of various operation conditions.

Product coode	212775
Mini bed attachment	GF200
Processing capacity	50 to 300g (It differs depending on whether the unit is of the batch type or specific samples used.)
Flow layer chamber capacity	3L
Spray nozzle	Dual fluid nozzle: 1A standard
Stirring blades	Integrated inside the flow layer chamber
Filter	Polyester (Carbon fiber mixed PTFE membrane laminate)
Filter cleaning mechanism	Pulse jet system
Glass parts	Ultra hard glass
Weight	Approx. 13 kg



Specifications

Product code	212778			
Model	GB210-B (GB210(basic unit)+GF200)			
1110 1101				
Temp. adjusting range	40 to 220°C (inlet temperature), 0 to 60°C (Outlet temperature)			
Temperature adjusting accuracy	Inlet temperature ± 1°C			
Spraying system	Two-way nozzle, Nozzle No. 1A as standard			
Drying air amount adjusting range	0 to 0.7m³/min			
Spray air pressure adjusting range	0 to 0.3MPa			
Liquid sending pump flow rate range	0 to 26mL/min			
External output	Inlet temperature, Outlet temperature, Temperature outlet (4-20 mA)			
Automatic lift	Moving up/down of glass chamber automatic lift			
Temperature adjusting device	PID digital temperature adjusting device			
Touch panel	Blower, Heater, Liquid sending pump, Pulse jet switch, Error display			
Control select switch	Inlet temperature, output temperature control switch (Outlet temp. control is conditional)			
Temperature sensor	K-thermocouple			
Heater	2.0 kW (at 200V) to 2.88 kW (at 240V)			
Liquid sending pump	Fixed amount peristaltic pump			
Spraying air pump	Spraying air compressor (Sold separately) is used			
Service outlet	For stirrer: AC100V, Max. 2A			
Suction blower	Bypass blower, Brushless DC motor			
Filter	Suction filter, Exhaust filter			
Spray nozzle cooling mechanism	Connector: Nipple×2, O.D.: ø10.5mm			
Spray air connection diameter	Nipple diameter: ø7mm			
Exhaust connecting diameter	ø50mm			
Safety device	Inlet/Outlet temperature overheat, Sample feed reverse rotation mechanism, Over current electric leakage breaker, Nozzle connection error			
External dimensions*1	W760×D420×H1,350 mm			
Weight	Approx. 110 kg			
Power supply (50/60Hz) rated current	AC220V 17A, AC240V 18A, Switching of terminals necessary			
Accessories	Silicon tube (with a stopper)×3, Tiron tube (with a stopper)×2, Exhaust duct (with one hose band)×1, Outlet temperature sensor, Spray air tube, Sample box, Static electricity removal earth, Teflon braided hose 5m (with two hose bands), Container table			
Necessary utility	28L/min. air volume and 0~294kPa(3kg/cm²) compressed air is required			
The length of the power gord is shout 2m outside the unit				

The length of the power cord is about 2m outside the unit.

- *1 External dimensions do not include projections.
- The air compressor used in this system must have a pressure regulator with air flow of 20 L/min or more and discharge pressure of 0 to 294kPa (3kg/cm²).
- Please note that this equipment is not explosion-proof for use with flammable or explosive substances.

Yamato Scientific Co., Ltd.

GB210-B

www.yamato-scientific.com
2003

Example of implementation

Sample Bind		Binder	inder		Test conditions				Results		
Name	Weight (min)	Name	Density (%)	Spray amount (min)	Inlet temp. (°C)	Liquid sending rate (g/min)					12~115 mesh recovery rate(%)
Silicon	200	PVA	5.0	77	125	15	59 (0.6)	4	27	339	58
Oxidized iron	160	PVA	2.5	50	120	15	98 (1.0)	4	21	205	62
Ceramics	200	PVA	3.0	106	120	15	78 (0.8)	3	22	404	82
Alumina	160	PVA	3.0	60	110	15	59 (0.6)	4	22	311	88
Silica	150	CMC	1.0	100	120	15	78 (0.8)	4	22	306	60
Lactose	200	Sorbitol	70.0	10	100	14	98 (1.0)	4	25	390	80
Black tea essence	250	Guar gum	0.5	24	85	6	59 (0.6)	10	28	333	77
Grease containing powder	200	Glucose	30.0	11	85	4	59 (0.6)	7	22	236	82

^{*}The average granule diameter is a geometric average.

Applications



Granulation, drying, mixing of powder Applications: Medicines, food, catalyst, die, detergent, ceramics, etc.

The unit accepts sample weight as less as 50 to 300g and is suitable for experiments of expensive samples or those of a labora-tory level.

Handling




Use of the one touch removal system has made removal or cleaning of the drying chamber, cyclone or the product container much easier.

Spraying Nozzle



Two-way nozzle system


The tip of the nozzle comprises of a nozzle for liquid and a nozzle for gas.



Model	Nozzle No.	Size (µm)	Particle size	
1A	(F) 1650	A 406 B 1270	1~40µm	
(Standard)	(A) 64	C 1626		
1	(F) 2050	A 508 B 1270	5~40µm	
-	(A) 64	C 1626		
2A	(F) 2050	A 508 B 1270	5~50µm	
	(A) 70	C 1778		
2	(F) 2850	A 711 B 1270	10~40µm	
	(A) 70	C 1778		
3	(F) 2850	A 711 B 1270	10~50µm	
	(A) 64	C 1626		

Particle sizes may vary on samples used and parameter settings.

Diagram



No.	Part name
1	Heater
2	Micro porous plate
3	Flow layer chamber
4	Filter chamber
(5)	Nozzle
6	Filter
7	Liquid sending pump
8	Blower
9	Interim pipe
10	Solenoid valve

1	No.	Part name						
	11)	3-way solenoid valve						
	12)	Needle valve						
	Pressure meter							
	14 Stirring blades							
	15)	Stirring motor						
	16 Inlet temperature sensor							
	17)	Outlet temperature sensor						
	18)	Blind (service port)						
	19	Suction port, suction filter						
	20	Nozzle cooling connection						

Optional items

	•	Desil of
No.	Product name	Product code
1	Safety cover	212787
2	Viton packing for cyclone inlet/outlet (1 set of 2 types)	212781
3	Teflon packing for cyclone inlet/outlet (1 set of 2 types)	212782
4	Regulator	212789
(5)	Supply air filter box (for 0.3 µm collection)	212791
6	*Inlet/outlet temperature recorder (3-dot)	212792

^{*}Please specify when ordering main unit